



ELECTRIC TITANIUM BED

3 MOTOR

MB105




EASY SIDE GUARD MECHANISM

Fitted with self locking mechanism to ease full adjusted side guard panel.



MATTRESS PLATFORM

Made of strong and lightweight material. Easy to clean.



AUTO REGRESS

The Auto Regress System reduce pressure in patient's pelvic area when the back-raise is raised.




INDIVIDUAL LOCK

Four 5" single wheel castors with individual locking system.

SPECIFICATIONS

Construction	Combination of high quality steel section, and sheet construction utilizing the latest manufacturing techniques throughout.
Head & Foot End Panel	Polypropylene panels (Blow molding process).
Side Guard	Polypropylene side guards, independent folding and interlocking system.
Positions	Back-raise inclination 0 - 70°. Knee-raise inclination 0 - 24°. Mattress platform elevation 400 (Lo) - 730 (Hi)mm CPR - Electric single touch button. Auto contour.
Castor	5" Single wheel castors, two with individual lock.
I.V. Pole	Stainless steel sections - adjustable. Can be located in four different positions.
Accessories	Accessories Railing, Backrest angle indicator.
Coating	Epoxy powder coating.
Patient Max. Load	185 kg.
Safe Working Load	250 kg.
Net Weight (approx.)	140 kg.
Dimension (approx.)	(L) 2220 x (W) 1040 x (H) 780 - 1115 mm (overall). (L) 1975 x (W) 900 x (H) 400 - 730 mm (on mattress platform).

DIMENSIONS (mm)



Product images are for illustrative purposes only and color may differ from the actual product.
Specifications subject to change without prior notice



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Eff. Date 19-01-2023



Marketing Office :
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com
www.mak-techno.com