


ELECTRIC LOVINA BED

3 MOTORS

MB445



STYLE AND SAFETY

The modern and full length side guard ensures patient safety. With its 370mm height from the mattress platform, the side guard will reduce the risk of falling.

SIDE GUARD

Lovina side guards have been specially designed in compliance with IEC 60601-2-52 to prevent entrapment risk.

AUTO-REGRESSION SYSTEM

The Auto Regress System reduce pressure in patient's pelvic area when the back-raise is raised.



MATTRESS PLATFORM

Made of strong and lightweight material. Easy to clean.

SPECIFICATIONS

Construction	Combination of high quality steel section, and sheet construction utilizing the latest manufacturing techniques throughout.
Head & Foot End Panel	Polypropylene panels (blow molding process).
Side Guard	Polypropylene side guards, independent folding and interlocking system.
Positions	Back-raise inclination 0 - 70°. Knee-raise inclination 0 - 24°. Mattress platform elevation 370(Lo) - 700(Hi) mm.
Castor	Four 5" single wheel castor, two with brakes.
Power Consumption	100-240V AC, 50/60Hz, 350W.
I.V. Pole	Stainless steel sections - adjustable. Can be located in four different positions.
Accessories	Inclinometer back-raise and Accessories rail.
Coating	Epoxy powder coating.
Max. Patient Weight	185 kg.
Save Working Load	250 kg.
Net Weight (approx.)	125 kg.
Dimension (approx.)	(L) 2200 x (W) 1060 x (H) 755 - 1085 mm (overall) (L) 1975 x (W) 900 x (H) 370 - 700 mm (on mattress platform).

DIMENSIONS (mm) ±



Product images are for illustrative purposes only and color may differ from the actual product
Specifications subject to change without prior notice



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Eff. Date 21-11-2023



Marketing Office :
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com
www.mak-techno.com