


# ELECTRIC TITANIUM BED


## 4 MOTORS

# MB106




### SIDE GUARD

Equipped with self locking mechanism, the side guards provide protection to patients.



### BED EXTENSION

Bed extension (up to 20cm) accommodates different patient sizes.



### AUTO REGRESSION

Auto regress provides maximum comfort for patient and reduce the risk of pressure ulcer.




### CPR QUICK RELEASE

For fast response in the event of CPR, this bed comes with manual and electric quick release.

## SPECIFICATIONS

<b>Construction</b>	Made with high quality steel section, and sheet construction utilizing the latest manufacturing techniques throughout.
<b>Head &amp; Foot End Panel</b>	Polypropylene panels (Blow molding process).
<b>Positions</b>	Back-raise inclination 0 - 70°. With Auto regress, approx. 140mm Knee-raise inclination 0 - 24°. Mattress platform elevation 400 (Lo) - 740 (Hi)mm, with bed extension. Trendelenburg - anti trendelenburg -14° - 14 ° CPR - Electric single touch button. Auto contour.
<b>Castor</b>	5" Single wheel castors, two with individual lock.
<b>I.V. Pole</b>	Stainless steel sections - adjustable. Can be located in four different locations.
<b>Accessories</b>	- Accessories Railing - Angle Indicator - X-ray Holder
<b>Optional</b>	Back-raise mattress board for x-ray function.
<b>Coating</b>	Epoxy powder coating.
<b>Max. Patient Weight</b>	185 kg.
<b>Safe Working Load</b>	250 kg.
<b>Net Weight (approx.)</b>	144 kg.
<b>Dimension (approx.)</b>	(L) 2220 x (W) 1070 x (H) 780 - 1120 mm (overall) (L) 1975 x (W) 900 x (H) 400 - 740 mm (on mattress platform)

## DIMENSIONS (mm)



Product images are for illustrative purposes only and color may differ from the actual product  
Specifications subject to change without prior notice

Eff. Date 19-01-2023

**EC REP** MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany



**Marketing Office :**  
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940  
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com  
www.mak-techno.com