PEDIATRIC BED

MBC01







CASTOR Giving smooth bed moving fitted 5" single wheel castor, 2 with brakes system

SPESIFICATIONS

Construction



PANEL

Fully steel sections, provide excellent looks and safety for patient during treatments.

Combination of high quality steel section, and sheet



SIDE GUARD RAIL

Full steel section side guard rail ensures provide excellent looks and safety for children during treatments. With sliding system enable full adjusted its side guard rail



COLLAPSIBLE CRANK Equipped excellent crank driven enable full adjusted product solution

DIMENSIONS (mm) ±



construction utilizing the latest manufacturing techniques throughout. **Mattress Platform** Manufactured from phosphate coated perforated steel sheet with mattress locators to sides and ends. Panel Head & Foot End High quality steel sections panel. Side Guard Sliding, steel sections side guard. Position Back-raise inclination 0° - 70°. - Manually by cranks Castor 5" swivel lock castors, two with brakes. I.V. Pole Stainless steel sections - adjustable. Can be located in four different positions. Finishing Epoxy powder coating. Max. Load 120 kg. 57 kg. Nett Weight (approx.) (L)1770 x (W)815 x (H)1390 mm (overall). **Dimensions** (approx.) (L)1700 x (W)700 x (H)690 mm (on mattress platform).

Product images are for illustrative purposes only and color may differ from the actual product. Specifications subject to change without prior notice

Eff. Date 19-01-2023



EC REP MDSS GmbH Schiffgraben 41 30175 Hannover, Germany



Marketing Office :

Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940 Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com www.mak-techno.com