



01

■ ELECTRIC DELIVERY BED 34191

Eff. Date: 19-03-2023



HI-LO POSITION

Height of mattress base 500 - 840 mm



Protection

TRENDELENBURG POSITION Trendelenburg 13°



BACKREST POSITION Backrest 0° to 70°

DIMENSION (mm) ±









CALF SUPPORT



PATIENT HANDLE



I.V POLE

Features DELIVERY BED



CPR HANDLE

Handle CPR is used for increase the effectiveness of CPR actions by providing support appropriate for compression optimal chest.



GADGET POCKET STORAGE

Additional storage on the side side guard that can used for saving patient's personal items.



MATTRESS EXTENSION

Setting with sliding mechanism and automatic locking as flexibility and effectiveness of the position for mothers during the birthing process.



CENTRAL LOCK CASTOR

Central lock castor provides ease of use with the ability to lock everything wheels simultaneously, improve stability and mobility product quickly.



BUMPER

The bumper protects the product from damage when it occurs light impact



HAND REMOTE

Complete comfort and control over the patient's sleeping position with innovative hand remote control.



BASKET

Container with a capacity of 9 liters which is safe and hygienic for place fluids during the process giving birth, its use can be removable.



Mattress with a width of 900 mm which is enough for comfort Mother and baby. Made from foam with quality synthetic leather and easy to clean.



SIDE GUARD

2 side guards to prevent damage the patient fell. Easy to open and close, just 1 hand.

EXCLUSIVE MATTRESS



PT MEGA ANDALAN KALASAN

MARKETING OFFICE

Gran Rubina Business Park Lt. 22C, Komp. Rasuna Epicentrum, Jl. HR. Rasuna Said, Jakarta, Indonesia 12940

t. (+62) 21 - 837 00 555 f. (+62) 21 - 837 00 335 marketing@mak-techno.com

International Contact: export@mak-techno.com

FACTORY

Jl. Tanjung Tirto 34, Tirtomartani, Kalasan, Km 13, Yogyakarta Indonesia 55571

t. (+62) 274 - 497 069 f. (+62) 274 - 496 226

www.mak-techno.com