


MOLUCCAS BED

2 CRANK

MB512



SMALL WIDTH
Suitable for narrow space.



PANEL
The head and foot end panel are easy to clean, lightweight and designed with self-locking system for easy assembly and disassembly.



CASTOR
5" castors with brakes.




COLLAPSIBLE SIDE GUARD
Side guard locking system with one touch locking mechanism to ensure patient safety.

SPECIFICATIONS

| | |
|----------------------------------|---|
| Construction | Manufactured with high quality steel and sheet construction. |
| Mattress Platform | Made from perforated steel sheet with mattress stopper at all sides. |
| Head & Foot End Panel | Polypropylene panels. |
| Side Guard | Collapsible, one-touch locking aluminum side guard. |
| Position | - Back-raise inclination : 0 - 70°. - Knee-raise inclination : 0 - 25°. Collapsible crank driven. |
| Castor | Four 5" single wheel castor, two with brakes. |
| IV. Pole | Stainless steel - adjustable. Can be located in six different positions. |
| Finishing | Epoxy powder coating. |
| Max. Load | 170 kg. |
| Safe Working Load | 200 kg. |
| Net. Weight (approx.) | 77 kg. |
| Dimension (approx.) | (L) 2155 x (W) 860 x (H) 835 mm (overall). (L) 1985 x (W) 770 x (H) 500 mm (on mattress platform). |

DIMENSIONS (mm) ±



Product images are for illustrative purposes only and color may differ from the actual product.
Specifications subject to change without prior notice



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Eff. Date 18-01-2023



Marketing Office :
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com
www.mak-techno.com