

# LA NINA BED 2 CRANKS

73062



## CRANK

Back-raise and Knee-raise positions is adjustable by foldable crank.



## PANEL

The head and foot end panel are easy to clean and designed with self-locking system for easy assembly and disassembly.



## MATTRESS PLATFORM

Made of strong and lightweight material. Easy to clean.




## COLLAPSIBLE SIDE GUARD

Side guard locking system with One Touch Locking Mechanism to ensure patient safety.

## SPECIFICATIONS

<b>Construction</b>	High quality steel, and sheet construction utilizing the latest manufacturing techniques throughout.
<b>Mattress Platform</b>	Made of polypropylene plastic with high quality steel frame.
<b>Head &amp; Foot End Panel</b>	Detachable ABS panel with steel frame.
<b>Side Guard</b>	Collapsible, one-touch locking aluminum side guard.
<b>Position</b>	Back-raise inclination 0-70°. Knee-rest inclination 0-24°. - (collapsible crank driven).
<b>Castor</b>	Four 5" single wheel castor, two with brakes.
<b>I.V. Pole</b>	Stainless steel - adjustable. Can be located in four different positions.
<b>Coating</b>	Epoxy powder coating finish.
<b>Max. Load</b>	170 kg.
<b>Net Weight (Approx.)</b>	75 kg.
<b>Dimension (Approx.)</b>	(L)2160 x (W)990 x (H)885 mm (overall). (L)1920 x (W)840 x (H)550 mm (on mattress platform)
<b>Packaging</b>	Knock-down system.

## DIMENSIONS (mm) ±



Product images are for illustrative purposes only and color may differ from the actual product.  
Specifications subject to change without prior notice

Eff. Date 07-05-2024



**Marketing Office :**  
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940  
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com  
www.mak-techno.com