


# BORNEO Jr. MANUAL BED


## 2 CRANK

# MB212



### CRANK

Backrest & kneerest position is adjustable by foldable crank



### PANEL

Head & foot end panel are easy to clean and designed with self-locking system for easy assembly & dissambely



### CASTOR

5" castor with brakes for easy mobility




### SIDE GUARD

Collapsible side guard with one touch locking system mechanism to ensure patient safety

## SPECIFICATIONS

<b>Constructions</b>	High quality steel and sheet construction utilizing the latest manufacturing techniques throughout.
<b>Mattress Platform</b>	Polypropylene.
<b>Head &amp; Foot Panel</b>	Detachable polypropylene.
<b>Side Guard</b>	Collapsible, one touch locking aluminium side guard.
<b>Positions</b>	- Backrest inclination 0 - 70°. - Kneerest inclination 0 - 25°. Collapsible crank driven.
<b>Castor</b>	Four 5" single wheel castor, two with brakes.
<b>I.V. Pole</b>	Stainless steel frame, adjustable. Can be located in four different positions.
<b>Finishing</b>	Epoxy powder coating.
<b>Net Weight</b>	± 71 kg.
<b>Max. Load</b>	170 kg.
<b>Safe Working Load</b>	200 kg.
<b>Dimension</b>	± (L) 2130 x (W) 1000 x (H) 835 mm (overall). ± (L) 1920 x (W) 840 x (H) 500 mm (on mattress platform).

### DIMENSIONS (mm) ±



Product images are for illustrative purposes only and color may differ from the actual product.  
Specifications subject to change without prior notice



MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany



Eff. Date 15-02-2024



**Marketing Office :**  
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940  
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com  
www.mak-techno.com