

MANUAL CELEBES BED

3 CRANKS (SIDE BLOW)

74003



STYLE & SAFETY

The modern and full length side guard ensures patient safety. With its 370mm height from the mattress platform, the side guard will reduce the risk of falling.



SIDE GUARD

Celebes side guards have been specially designed in compliance with IEC 60601-2-52 to prevent entrapment risk.



AUTO-REGRESS

The Auto Regress System reduce pressure in patient's pelvic area when the back-raise is raised.



MATTRESS PLATFORM


Made of strong and lightweight material. Easy to clean.

SPECIFICATIONS

Construction	Combination of high quality steel section, and sheet construction utilizing the latest manufacturing techniques throughout.
Head & Foot End Panel	Polypropylene panels (Blow molding process).
Side Guard	Polypropylene side guards, independent folding and interlocking system, with inclination indicator.
Positions	Back-raise inclination 0 - 70°. Knee-raise inclination 0 - 25°. Mattress platform elevation 370 (Lo) - 680 (Hi) mm. - Manually by cranks. (torque release feature).
Castor	5" single wheel castor, 2 with brakes.
I.V. Pole	Stainless steel sections - adjustable. Can be located in four different positions.
Accessories	Backrest angle indicator and accessories railing.
Coating	Epoxy powder coating.
Max. Patient Weight	170 kg.
Safe Working Load	200 kg.
Net Weight (Approx.)	123 kg.
Dimension (Approx.)	(L) 2200 x (W) 1060 x (H) 755-1065 mm (overall). (L) 1975 x (W) 900 x (H) 370-680 mm (on mattress).

Specifications subject to change without prior notice

DIMENSIONS (mm)



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Eff. Date 10-05-2024



Marketing Office :
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com
www.mak-techno.com