


BORNEO Jr. MANUAL BED


3 CRANK CLC

MB214



CRANK

Backrest, kneerest, & hi-lo position is adjustable by foldable crank



PANEL

Head & foot end panel are easy to clean and designed with self-locking system for easy assembly & dissambely



CASTOR

5" castor with brakes for easy mobility




SIDE GUARD

Collapsible side guard with one touch locking system mechanism to ensure patient safety

SPECIFICATIONS

Constructions	High quality steel and sheet construction utilizing the latest manufacturing techniques throughout.
Mattress Platform	Polypropylene.
Head & Foot Panel	Detachable polypropylene.
Side Guard	Collapsible, one touch locking aluminium side guard.
Positions	<ul style="list-style-type: none"> - Backrest inclination 0 - 70°. - Kneerest inclination 0 - 25°. - Mattress platform elevation 450 (Lo) - 650 (Hi) mm. Collapsible crank driven.
Castor	Four 5" double wheel castor, with central lock mechanism.
I.V. Pole	Stainless steel frame, adjustable. Can be located in four different positions.
Finishing	Epoxy powder coating.
Net Weight	± 95 kg.
Max. Load	170 kg.
Safe Working Load	200 kg.
Dimensions	± (L) 2130 x (W) 1000 x (H) 785 - 985 mm (overall). ± (L) 1920 x (W) 840 x (H) 450 - 650 mm (on mattress platform).

DIMENSIONS (mm) ±



Product images are for illustrative purposes only and color may differ from the actual product.

Specifications subject to change without prior notice



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Eff. Date 15-02-2024



Marketing Office :
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com
www.mak-techno.com