


# LA NINA BED

## 3 CRANKS CLC

73064



### CRANK

Equipped with excellent crank driven to adjust mattress position.



### PANEL

The head and foot end panel are easy to clean and designed with self-locking system for easy attachment.



### CENTRAL LOCK

Four 5" double wheel castors with central locking system.



### COLLAPSIBLE SIDE GUARD


Side guard locking system with One Touch Locking Mechanism to ensure patient safety.

## SPECIFICATIONS

<b>Construction</b>	High quality steel and sheet construction utilizing the latest manufacturing techniques throughout.
<b>Mattress Platform</b>	Made of perforated polypropylene lining with high quality steel sheet frame.
<b>Head &amp; Foot End Panel</b>	Detachable ABS panel with steel frame.
<b>Side Guard</b>	Collapsible, one-touch locking aluminum side guard.
<b>Position</b>	Back-raise inclination 0-70° Knee-rest inclination 0-24°. Mattress platform elevation 450 (Lo) - 650 (Hi) mm - (collapsible crank driven).
<b>Castor</b>	Four 5" double wheel castors - central locking system.
<b>I.V. Pole</b>	Stainless steel - adjustable. Can be located in four different positions.
<b>Coating</b>	Epoxy powder coating finish.
<b>Max. Load</b>	170 kg.
<b>Safe Working Load</b>	220 kg.
<b>Net Weight (Approx.)</b>	97 kg.
<b>Dimension (Approx.)</b>	(L)2160 x (W)990 x (H)785 - 985 mm (overall). (L)1920 x (W)840 x (H)450 - 650mm (on mattress platform).

Specifications subject to change without prior notice

## DIMENSIONS (mm)



Eff. Date 04-05-2024

**EC REP** MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany



**Marketing Office :**  
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940  
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com  
www.mak-techno.com