



PREMIUM COUCH

CHANGING TABLE

34118



COLOR OPTION



SIDE GUARD

Side mattress covered with vinyl / synthetic leather.



FOOT SWITCH

Optional control for height adjustment.



LOCKING SYSTEM

Fitted with retractable casters and central lock mechanism to ease product's mobility.




FOOT STABILITY

Adjustable foot mechanism for stability on uneven surfaces.

SPECIFICATIONS	
Constructions	High quality steel sections, and sheet construction utilizing the latest manufacturing techniques throughout.
Upholstery	Manufactured in line using high quality foam, lining with vinyl/imitation leather.
HeightAdjustment	Using high-grade electric actuator motor Hi-Lo upholstery platform elevation 500 - 950 mm.
Side Guard	Upholstery mattress side guard.
Features	- Foot mechanism stability. - Hand control switch.
Coating	Epoxy powder coating finish.
Max. Load	200 kg.
Net Weight	± 67 kg.
Dimension	± (L)1970 x (W)610 x (H)500 - 950 mm.

DIMENSION (mm)



Ⓞ Castor for product movement only, do not move couch with patient on top.

Product images are for illustrative purposes only and color may differ from the actual product.
Specifications subject to change without prior notice.

Eff. Date 15-06-2023

EC REP MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Marketing Office :
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com
www.mak-techno.com