## **EXAMINATION TABLE**

## ME101





MATTRESS Wide mattress covered with synthetic leather, easy to clean



**BACK-REST** 

Adjustable padded backrest to provide patient comfort in treatment procedures



## FRAME

Made of high quality steel, strong and rigid



## FOOT PAD Prevent scratches or damaged on foot frame

SPESIFICATIONS	
Construction	Combination of high quality steel section, and sheet construction utilizing the latest manufacturing techniques throughout.
Upholstery	Manufactured in line with medical devices directive using high quality foam, laminated with synthetic leather.
Back Raise	Inclination 0° - 60° - telescopic mechanism.
Max. Load	170 kg.
Nett Weight (approx.)	± 33 kg.
Dimensions (approx.)	± (L)1920 × (W)800 × (H)735 mm.
Packaging	Knock-down

Specifications subject to change without prior notice





Schiffgraben 41 30175 Hannover, Germany



Eff. Date 25-01-2021

Marketing Office : Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940 Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com www.mak-techno.com

MDSS GmbH

EC

REP